

Pelletizer



Pelletizer Take-Off

The Pelletizer is uniquely designed for laboratory or pilot plant extrusion runs. It is offered as either a table-top unit or mounted on a mobile cart. The unit performs continuous, direct pelletizing of thermoplastic materials - with multiple extruded profiles, 1-5mm in diameter. Supplied with manual air cooling device.

Features:

- 1/17 HP DC Drive Motor for the Feed Roll
- 3/4 HP Motor for the Drive rotor
- Manual Air Cooling Device

It starts with a Brabender®...

...and ends with your success!!