


Order now!  
By phone or e-mail

## Calibration Kit:

- Farinograph®
- Extensograph®
- Amylograph®
- Viscograph®
- Micro Visco-Amylo-Graph®



Dedicated Service Line:  
**201-343-8425** Ext.99

## Calibration Kit for:

- Farinograph®
- Extensograph®
- Amylograph®
- Viscograph®
- Micro Visco-Amylo-Graph®

## Brabender® Calibration Kit

Peace of mind for you  
and your customers



**C.W. Brabender®**

50 East Wesley Street  
South Hackensack, NJ.  
07607

Tel: 201-343-8425  
servicedept@cwbrabender.com  
www.cwbrabender.com



35002-DE / Brabender CS / Calibration-Kit / 03 / 2016



... where quality is measured.



## Put your Brabender® instrument to the test!

Utilizing Brabender® reference material will enable you to benchmark actual and target values.

For reference checks, we supply our new Calibration Kit. The kit consists of our reference material tested against our master instrument in Duisburg, a master curve and a standard test method.



SAMPLE

Prepare your sample as described.



ANALYZE

Carry out your measurement.



COMPARE

Compare your readings with the reference curve and stated tolerances.



CHECK

If your values are within the permitted range, then your instrument works properly.

If your readings are unsatisfactory, please do not hesitate to contact us. We will check your instrument and find a solution.

### Fast and inexpensive

Check your instrument yourself using Brabender® reference material.

### Save money by taking pre-emptive action

Avoid expensive complaints, unnecessary waste or damage to your reputation by conducting regular checks.

### Certified safety

Show your customers that you take quality seriously: participate in our reference material program and regularly check your instruments.

